

DSM 9000 SMS Alarm System

(GSM Auto Dialer for
Temperature Alarm Alert
High & Low, Power Failure)


Laboratory Freezer & Refrigerator Connection

Features Include:

- ❖ Able to send SMS message when power failure or high/low temperature alarm alert is triggered.
- ❖ Max. link-up 7 units freezer(s)/refrigerator(s)
- ❖ Use Sim card device
- ❖ With external Audio & Visual alarm
- ❖ Remote alarm contact
- ❖ DC backup battery last for 72-hour
- ❖ Able to send max. 8 telephone numbers via SMS



Freezer Data Logging & Monitoring System v1.0



ETS BIO FREEZE SDN BHD (717189-H)

Add: No. 32G & 32-1, Dataran Mestika, Jalan Mawar 17, Taman Mestika,
56100 Kuala Lumpur, Malaysia.

Tel: +603 - 4292 2722 / 4288 2722

Fax: +603 - 4293 2722

Email: etsbiofr@gmail.com / sales_ets@hotmail.com

Website: www.etsbiofreeze.com



NOTE: In the event slight difference occurs between pictures and actual products, please refer to the actual products. We reserve the rights of final interpretation to this catalog and please contact us for future information you might need.