



DSM 9000 SMS Alarm System

(GSM Auto Dialer for
Temperature Alarm Alert
High & Low, Power Failure)


Laboratory Freezer & Refrigerator Connection

Features Include:

- ❖ Able to send SMS message when power failure or high/low temperature alarm alert is triggered.
- ❖ Max. link-up 7 units freezer(s)/refrigerator(s)
- ❖ Use Sim card device
- ❖ With external Audio & Visual alarm
- ❖ Remote alarm contact
- ❖ DC backup battery last for 72-hour
- ❖ Able to send max. 8 telephone numbers via SMS



Freezer Data Logging & Monitoring System v1.0



ETS BIO FREEZE SDN BHD (717189-H)

Add: No. 32G & 32-1, Dataran Mestika, Jalan Mawar 17, Taman Mestika, 56100 Kuala Lumpur, Malaysia.

Tel: +603 - 4292 2722 / 4288 2722

Fax: +603 - 4293 2722

Email: etsbiofr@gmail.com / sales_ets@hotmail.com

Website: www.etsbiofreeze.com

