



DSM 9000 SMS Alarm System

(GSM Auto Dialer for Temperature Alarm Alert, High & Low, Power Failure)


Laboratory Freezer & Refrigerator Connection

Features:

- ✓ Able to send SMS messages when power failure or high/low temperature alarm is triggered
- ✓ Max. link up 5 units freezer(s)/refrigerator(s)
- ✓ Use Sim Card device
- ✓ With external Audio & Visual alarm
- ✓ Remote alarm contact
- ✓ DC backup battery last for 72 hour
- ✓ Able to send max 5 telephone numbers via SMS



Freezer Data Logging & Monitoring System v1.0



ETS BIO FREEZE SDN BHD (717189-H)

Add: No. 32G & 32-1, Dataran Mestika, Jalan Mawar 17, Taman Mestika, 56100 Kuala Lumpur, Malaysia.

Tel: +603 - 4292 2722 / 4288 2722 / +6011-1077 4983

Fax: +603 - 4293 2722

Email: info@etsbiofreeze.com / sales_ets@hotmail.com

Website: www.etsbiofreeze.com



Contact Us

NOTE: In the event slight difference occurs between pictures and actual products, please refer to the actual products. We reserve the rights of final interpretation to this catalog and please contact us for future information you might need.